



Proseers

Portable Water Purification System



QINGDAO INNOVA BIO-MEDITECH CO., LTD





About Us

Innova Bio-Meditech is one of the leading global providers of laboratory and medical devices. Firmly committed to our mission of "sharing innovative bio-meditech solutions with the world" , we are dedicated to innovation in the fields of Biology Project, Life Science, Pharmacy Industry and Medical Treatment.

Innova Bio-Meditech possesses a sound distribution and service network with business partners in North and Latin America, Europe, Africa and Asia-Pacific etc, We have built up a well established R&D, manufacture network with 3 centers in Beijing, Qingdao, and shanghai. Inspired by the needs of our customers, we adopt advanced technologies and transform them into accessible innovation. This means constant effort and research, in order to more fully understand the developments of the market, produce constantly upgraded product ranges by adding new products year after year.


The passion for science

Content

Water purification phase	01
Applications of RO water	01
Applications of Ultra pure water	02
Water system selection guide	03
Overview	04
Smart and Unlimited	06
Unlimited expansion	06
Flexible access pure water	07
Excellent Water Quality	08
Purification cartridge	08
RO system	09
Dual wavelength UV lamp	10
TOC Monitoring	11
Optional point-of-use filter	11
Convenient Water Collect	12
Consumable identification	13
Advantages & Features	14
Process & Specifications	15
Dimensions & Operating Weights	17
Accessories	19
	20



Water purification phase



Applications of RO water

Reverse osmosis water (RO water) has a wide range of uses that cover many areas, the following are its uses:



RO water is used to clean semiconductors and circuit boards to ensure that there is no ionic contamination, which is critical to the performance of electronic products.

RO water is used in the pharmaceutical industry for drug preparation and equipment cleaning in the manufacturing process to ensure the purity and safety of drugs.



Experimental water: In scientific research, especially in chemical and biological experiments, RO water provides the necessary high-purity water source to guarantee the accuracy of experimental results.
 Precision Instruments: High precision analytical instruments used in laboratories often require RO water to generate or dilute standard solutions to avoid background interference.

Hospitals use RO water to clean surgical instruments to ensure a sterile environment and reduce the risk of post-operative infections.



Applications of Ultra pure water

Ultrapure water has been rigorously treated to remove almost all impurities, following are some of its major applications in several fields:



- ⦿ **Wafer cleaning:** In the semiconductor manufacturing process, ultra pure water is used to clean silicon wafers to ensure that their surfaces are free of particles and ionic contamination, thus ensuring high chip quality.
- ⦿ **Chemical solution dilution:** ultra pure water is used in the dilution process of photoresist and other chemicals to avoid the impact of impurities on the process.
- Equipment cooling: Many high-precision equipment needs to be cooled with ultra pure water to prevent corrosion and scale buildup.

- ⦿ **Drug preparation:** Ultra pure water is used in the preparation of injections, eye drops and other sterile preparations to ensure the high purity and safety of drugs.
- ⦿ **Laboratory analysis:** In drug development, ultra pure water is used for various experiments and analysis, such as high-performance liquid chromatography (HPLC).



- ⦿ **Cell culture:** In the process of cell culture, ultra pure water is used to formulate culture medium to avoid cell contamination.
- ⦿ **Molecular biology experiments:** Such as PCR (polymerase chain reaction), DNA sequencing and other experiments using ultra pure water to ensure the accuracy of the experimental results.

- ⦿ **Dialysis water:** Ultra pure water is needed during hemodialysis to avoid patient infection and complications.
- ⦿ **Pathology testing:** Ultrapure water is used to prepare reagents and clean samples to ensure the accuracy of test results.



- ⦿ **Analytical instruments:** High-precision instruments such as Inductively Coupled Plasma Mass Spectrometer (ICP-MS), Atomic Absorption Spectrometer (AAS), etc. require the use of ultra pure water to avoid background interference.
- ⦿ **Sample preparation:** Ultrapure water is used to prepare high purity samples in various chemical and biological experiments.

- ⦿ **Circuit board cleaning:** In the printed circuit board (PCB) manufacturing process, ultra pure water is used for cleaning to prevent residues from affecting circuit performance.
- ⦿ **High pressure boiler feed water:** In some demanding industrial applications, ultrapure water is used as boiler feed water to prevent scaling and corrosion.



Water system selection guide

Model	 Nimble RU8	 Nimble U
Feed water source	Potable tap water	Type III water
Provide RO water	✓	✗
Provide ultra pure water	✓	✓
Pre-guard cartridge 	✓	✗
Reverse osmosis 	✓	✗
RO water tank 	✓	✗
Ultra-pure UV 	✓	✓
Ultra purification cartridge 	✓	✓

Overview

|| Main structure



Front view

5-inch touch screen, displays on-line operation status, water quality information, consumables status and alarm information.

6 types of specialized terminal filters for different laboratory requirements.

5L portable separated design water tank, can be moved with one hand.



Left view

All 6 consumables can be recognized by NFC, identified with a single swipe.

Independent reverse osmosis design, ion retention rate of $\geq 98\%$.

Independent purification cartridge, filled with electronic-grade resin, can remove trace ions.




Right view

185/254nm dual wavelength UV lamps, effectively remove organic and bacteria.

*The model shown is Nimble RU8

Overview

Water leakage detector



External water leakage detector



Built-in water leakage detector

Independent water leakage protection, dual-site detection inside and outside the main host, accurate identification of water contact points, distinguish between water source leakage and instrument internal leakage, to ensure the safety of the entire laboratory. Once water leakage detected, 2 seconds after will turns to standby mode automatically.

RO water tank



5L RO water tank is easy to lift to each laboratory room , can be configured with multiple.

The water tank liquid level is monitored by pressure level sensor, and the liquid level display is accurate and stable; The storage capacity can be adjusted from 5% to 100% online.


Smart and Unlimited

Unlimited expansion

Flexible to satisfy your need

Choose the configuration that best suits your space and user needs, and then easily expand more modules as your lab expansion.

You can carry RO water tank of Nimble RU8 to Nimble U, and get ultra pure water by Nimble U.



The above configuration consists of 1x Nimble RU8 system and 3x Nimble U system

Smart and Unlimited


Flexible access pure water



1 Produce RO water. The RO water is stored in the RO water tank.



2 Carry the RO water tank.



3 Move to your lab room.



4 Connect the RO water tank to Nimble U. Dispense freshly ultrapure water.



5 Bring the RO water tank back to Nimble RU8. Produce RO water again.

The new Nimble series water purification system is simple, flexible and convenient, thus greatly improving your experimental efficiency.

Smart and Unlimited

Flexible access pure water



RO tank adapts to all Nimble series, and can produce ultra-pure water flexibly.

You can easily access freshly purified ultrapure water at every lab bench, even remote benches without raw water.

Excellent Water Quality

Purification cartridge

Pre-guard cartridge B

Filled with an appropriate amount of silicon phosphorus crystals to effectively reduce the hardness of feed water, Filled with wire wound filter elements to retain powder and floc impurities.

Ultra purification cartridge B

Filled with electronic grade ion exchange resin, the ions in the water are controlled at ultra-trace levels. Advanced vertical flow purification method is adopted to ensure the service life and purification effect of purification cartridge, and reduce the use cost.



*The model shown is Nimble RU8


Excellent Water Quality

Reverse osmosis system

Independent reverse osmosis design, the retention rate of reverse osmosis membrane in non-destructive state is $\geq 98\%$.

Water pumps

Qualified water pumps, low noisy and stable working pressure, ensure the stability of reverse osmosis; With intelligent feedback adjustment function, automatically adjust the working state of the water pump according to changes in temperature and water quality, ensure the stability of the system water production.



Excellent Water Quality

Dual wavelength UV lamp

185nm wavelength can oxidize organic compounds. 254nm wavelength can cause microbial DNA damage and has disinfection effect.



Total organic carbon Monitoring

Total organic carbon (TOC) data will be shown on the screen, evaluated by a built-in TOC monitoring. TOC level of Ultra pure water is always less than 5ppb($\mu\text{g/L}$).



TOC Monitoring

Excellent Water Quality

Optional point-of-use filter

You can choose to install 6 kinds of special terminal filters of the same brand to ensure that the water quality meets the requirements of various laboratory applications.



Standard	MIC-Filter	0.22µm filter; for removal of bacteria and particle(>0.22µm)
Optional	Bac-Filter	0.22µm+0.45µm filter; for removal of bacteria and particle(>0.22µm)
	Pyr-Filter	For removal of bacteria, protease, endotoxin, nuclease
	VOC-Filter	For removal of volatile organic compound-related applications
	ULC-Filter	For removal of ultra-trace organic matter
	EDS-Filter	For removal of endocrine disruptors

Convenient Water Collect

Click button



Click the button to collect ultrapure water, convenient and quick

Quantitative water dispense



Quantitative water dispense. Click the flask icon, enter the amount in milliliters, save and press "dispenser" to collect water, the collecting can be terminated by clicking "cancel"

RO water tank



Pour the RO water from the tank through the hole in the tank lid

Consumable identification

The consumables of this equipment include two purification cartridges, RO module, dual wavelength UV lamp, air filter and terminal filter, all six of which can be identified using NFC (base on RFID).

After entering the consumables identification interface, contact the chip on the consumables with the identification point on the side of the equipment, and the consumables information can be automatically identified.



The following consumables are all NFC (base on RFID) identifiable :



Reverse osmosis membrane



Ultra-pure UV



Purification column



Air filter



Point-of-use filter

Advantages & Features

Smart design system



1+N design, support module expansion.



Water tank separate design.



Multi-stage purification process combination.



A variety of water collect operations.




Whole tank transfer, no environmental interference.



Intelligent water injection module.

Cloud monitoring



Simultaneous monitoring of multiple devices from one account.



Consumables and alarm information are uploaded and monitored in real time.



Water quality data real-time monitoring, regular upload .



II Ultra purification



Medical polymer material synthesis, no dissolution.



Controlled particle size and pore size distribution for targeted removal of organic matter.



Electronic polishing resin, large work exchange capacity.



Gel state, no heat source and endotoxin, avoid the secondary pollution of resin.



II Intelligent human-computer interaction



5 inch touch screen, three languages display.



Three-level authority management.



Quantitative water collect: 0.02-60L adjustable.



Water quality reporting function, real-time monitoring.




Automatical data store of water quality, instrument maintenance and others.




Process & Specifications

|| Nimble RU8




|| Nimble U



Solenoid valve
 Electrode
 Air filter
 Flowmeter
 Waste-water ratio

Three-way valve
 Pump
 Pressure sensor
 Check valve
 Effluent liquor



Reverse osmosis membrane



Purification column



Ultra-pure UV



Pure water tank

Main configuration	All-in-one system	Ultra pure water
	Nimble RU8	Nimble U
Pre-guard cartridge	Yes	No
Reverse osmosis	Yes	No
Water tank	Yes	No
Ultra-pure UV	Yes	Yes
Ultra purification cartridge	Yes	Yes

Model	All-in-one system	Ultra pure water
	Nimble RU8	Nimble U
Feed water requirements		
Source	Potable tap water	Type III water
Water inlet pressure	0.1-0.4MPa	0-80KPa
Water inlet conductivity	<2000uS/cm	<100 us/cm
Water inlet TOC	<1ppm	<50 ppb
Ambient temperature	4-45°C	4-45°C
PH	4-10	6-8
Type I ultra purified water		
Electrical resistivity	18.2 MΩ.cm@25°C	
TOC (Total Organic Carbon)	<5 ppb ^[1]	
Particles (>0.2µm particles)	<1 /mL	
Microorganism	<0.01 CFU/mL	
Pyrogen	<0.001 EU/mL	
Rnase	<0.4 pg/mL	
Dnase	<2 pg/mL	
Flow rate	1L/min	
RO water		
RO rejection	97-98%	N/A
Organic rejection	≥99%	
Particulate and bacterial rejection	> 99%	
Electrical requirement		
Input voltage	100-240V/50-60Hz/AC	100-240V/50-60Hz/AC
Rated power	120w	120w
Size information		
Dimensions (width × depth × height)	Main host: 220*510*445 mm Water tank: 135*260*400 mm	Main host: 220*510*445 mm
Main host weight	12 kg	9 kg

[1] Under proper operating conditions


Dimensions & Operating Weights



Nimble RU8



Water tank



Nimble U



Indoor use only.



Avoid direct sunlight.



Altitudes between 0 and 2000 meters.




Ambient temperature: 4~45°C.



Operating voltage: 100~240V/50~60/AC.



A clear drainage pipe 1.5m from the instrument.



Fixed, horizontal, and non-combustible mounting surface.



No heat source around



Away from strong magnetic fields.



Accessories

Purified parts



Reverse osmosis membrane



Ultra-pure UV



Purification column

Point-of-use filter



MIC-Filter



Bac-Filter



Pyr-Filter



VOC-Filter



ULC-Filter



EDS-Filter

Others



Air filter



Water leakage detection



Proseers



Qingdao Innova Bio-Meditech Co., Ltd.

Add.: No.176 Jufeng Road,266199, Qingdao, China

Tel.: +86 532 8789 0634

Email: info@innobiomed.com

Web: www.innovabiomed.com